



**Forum for  
Collaborative HIV Research**

# **LOD/LLOQ Discussion 7<sup>th</sup> HCV DRAG November 3<sup>rd</sup>, San Francisco**

Gabrielle Heilek Ph.D.  
Director Medical Affairs and Clinical Research Lead,  
Hepatitis  
Roche Molecular Diagnostics

on behalf of the RMD Medical/Scientific Affairs Team



# Definition of Terms

- **Target not detected (TND):**
  - No PCR amplification
  - Equivalent to **0 IU/ml** HCV RNA
  - Validated via Analytical Specificity studies
- **Limit of Detection (LOD)**
  - Lowest amount of analyte that can be detected 95% of the time as a positive PCR signal
  - Varies with matrix and genotype
  - **Not** reported on validated test software
- **Lower limit of Quantification (LLOQ):**
  - Lowest amount of analyte to be both linear and accurate
  - Validated in linearity and total analytical error studies

# Commercially Available HCV RNA PCR Tests

*LOD and LLOQ assay performance*

Assay	Linear Range	LLOQ	Overall LOD <sup>1</sup>
COBAS TaqMan HCV Test, v2.0 for use with the High Pure System	25 – 300,000,000 IU/mL	25 IU/mL	20 IU/mL
COBAS AmpliPrep/COBAS TaqMan HCV Test Version 1	43 – 69,000,000 IU/mL	43 IU/mL	CE-IVD 15 IU/mL US-IVD 18 IU/mL
VERSANT® HCV RNA 3.0 Assay (bDNA) <sup>2</sup>	615 - 7,690,000 IU/mL	615 IU/mL	1,000 IU/mL
Abbott RealTime HCV Test	12 IU – 8.00 Log <sub>10</sub> IU/mL	12 IU/mL	12 IU/mL
COBAS AmpliPrep/COBAS TaqMan HCV Test Version 2 <b><i>In Development</i></b>	<i>Targeted Range: 15 - 8.00 Log<sub>10</sub> IU/mL</i>	<i>LOD=LLOQ</i>	

<sup>1</sup> Performance for LOD is based on overall genotype data

<sup>2</sup> Assay LOD at 1,000 IU/mL (95% detection rate), LLOQ and range determined per package insert, 615 IU/mL is <95% detection rate

# High Pure TaqMan HCV v2 Test Results

*Validated software read out and lab report guidance*

Titer Result	Results Interpretation
Target not Detected	Report results as “ <b>HCV RNA not detected</b> ”
< 25 IU/ml	Calculated IU/ml is below the lower limit of quantification (LLOQ) of the assay. Report results as “ <b>HCV RNA detected, &lt;25IU/ml, HCV RNA is not quantifiable</b> ”
$\geq 25$ IU/ml and $\leq 3.0E+08$ IU/ml	A <b>quantifiable number</b> in the linear range of the assay
$> 3.0E+08$ IU/ml	Results are above the range of the assay. Report results as “ <b>&gt; 3.0E+08 IU/ml HCV RNA</b> ”

\* FDA approved COBAS® TaqMan® HCV Test, v2.0 for use with the High Pure System, package insert


Assay used during clinical studies for boceprevir and telaprevir.  
Gained FDA approval and global commercialization

# Conclusion and Recommendation

- Lab should report out only validated software prompts
  - No combination of TND and <LLOQ with fusion of terms
- LOD not suitable as a clinical parameter as it varies with country/regional approval, matrix and genotypes
- Clinical Trial analysis:
  - **TND** and **fixed cut-off** within linear range ( i.e. 25 or 50 IU/ml FDA call June 30th 2011) of several assays most suitable for SVR predictions

# Viral Load Decline on Therapy

*What does this all mean for drug development and clinical practice?*



# Backup

# High Pure Taqman HCV v2 Test Results

*The FDA approved LOD in the US-IVD Package Insert*

*plasma*

Genotype	Geometric Mean of Observed Results at the Lowest Concentration with > 95% Positivity Rate	Number of Replicates Tested	Number of Positive Results	Positivity Rate
1	15.1 IU/mL	66	66	100%
2	5.6 IU/mL	137	131	96%
3	15.3 IU/mL	63	63	100%
4	12.0 IU/mL	63	63	100%
5	3.7 IU/ml	63	60	95%
6	20.4 IU/mL	66	65	98%

LOD differs

- depending on *sample* used (*plasma, serum*)
- between *genotypes 1-6*

*serum*

Genotype	Geometric Mean of Observed Results at the Lowest Concentration with > 95% Positivity Rate	Number of Replicates Tested	Number of Positive Results	Positivity Rate
1	2.1 IU/mL	66	63	95%
2	10.9 IU/mL	126	124	98%
3	2.5 IU/mL	68	65	96%
4	7.5 IU/mL	71	71	100%
5	12.1 IU/mL	66	65	98%
6	7.9 IU/mL	64	62	97%

**Limit of Detection was assigned using highest titer with >95% positivity rate and the most conservative genotype: 20 IU/ml**



# Why LOD is not a General Usable Term

## *CE vs. US Registration Data Differences*

Genotype 1	CE Registration Data	US Registration Data
Plasma	9.3 IU/ml	15.1 IU/ml
Serum	8.8 IU/ml	2.1 IU/ml

N>148 per matrix and concentration

N>63 per matrix and concentration