




**THE FORUM**  
For Collaborative Research<sup>SM</sup>

# FibroNest: Single-fiber, High Content Digital Pathology and AI

Mathieu M. Petitjean, PhD

CEO, PharmaNest, INC, Princeton, USA



# About the FibroNest method

- ❑ **Same Slide(s) as Pathologist**
- ❑ **High content |Single Fiber / Nuclei**
- ❑ **AI > large quantitative & Relevant data-lakes**
- ❑ **Fully translational, not trained in existing paradigms**
- ❑ **1 Phase 3, 6 Phase 2 NASH, 3 Phase 2 in other Fibrotic conditions, 60+ preclinical studies**






**THE FORUM**  
For Collaborative Research™

# FibroNest: Fibrosis

# and Inflammation



- Normal Hepatocytes
- Steatosis Hepatocytes
- Inflammatory Cells
- Specialized Cells
- Q2-2022
- Glycogenotic Nuclei
- Acidophilic Bodies
- Megamitochondria
- Q2-2023
- Stellate Cells
- FibroBlasts





# Phenotypes of Fibrosis

Scar

Opportunity

Scar

Opportunity

Opportunity

Scar

Fiber Density (Faint-Mature)



Architecture Complexity



Macro Steatosis (Medium-large)





# Phenotypes of Fibrosis

Opportunity →

← Opportunity

Opportunity →


Fiber Density (Faint-Mature)

Architecture Complexity

Macro Steatosis (Medium-large)




# FibroNest quantifies Response to Treatment (Aramchol, NCT02279524)



Biopsy methodology	Post-BL Biopsy at <W48 weeks		Post-BL Biopsy at ≥ W48	
	N	%	N	%
All	28	100%	23	100%
Fibrosis Improvement (1 point or more) based on NASH CRN	7	25%	9	39%
Fibrosis Improvement (Paired reading ranked assessment) based on comparing individual patients slides	12	43%	14	61%
Subject Fibrosis Response (AI reading) using FibroNest's Phenotypic FCS (A responder is defined by an absolute reduction of > 0.3)	15	54%	23	100%
Subject Fibrosis Response (AI reading) using FibroNest's Phenotypic FCS (A responder is defined by a relative reduction of 25%)	6	21.4%	15	65.2%


FibroNest Parenchymal continuous Resolve paired-biopsy drug effect “inside Categorical Scores” and decouples anti steatosis from anti-fibrotic dynamics



## F4 Phenotypes


- Septal Region Morphometry
- Lobular Tissue Morphometry
- “Scar” Phenotypic Score to balance Tissue

## Injury and Tissue Regeneration



## F4 Sub phenotyping

- ASH-NASH
- LAENNEC (WIP)




## F4 Outcomes Prediction

- HCC progression
- Clinical Outcomes Prediction (WIP)




## Histological Parading Responders




Improved: 5 (29%)  
Remained stable: 8 (47%)  
Worsened: 4 (24%)

## Fibrosis Single Fiber Phenotyping



## POISE Study – OCA in PBC patients ( N=17)

## Continuous Scoring Response



25% Mean change from baseline  
51% responders




# PSC Opportunity

- ❑ **Histological Phenotype Fibrosis**
- ❑ **Continuous scores for**
  - **Fibrosis Severity**
  - **Progression risk**
  - **Phenotypes of treatment response**
- ❑ **Unsupervised ML. No need for annotations**
- ❑ **Same stained slides as pathologists**
- ❑ **N=20 for POC**



# Thank You !

Mathieu.Petitjean@pharmanest.com



PHARMANEST  
Princeton, USA

[www.fibronest.com](http://www.fibronest.com)