

Method comparison for evaluating new (improved) CD4 and viral load assays for laboratory service

Lesley Scott, PhD

Department of Molecular Medicine and
Haematology, University of the
Witwatersrand, Johannesburg,
South Africa




Why perform evaluations?

- **Justify** or **confirm** before routine implementation (result reporting)
- Sensitivity and specificity in local population (VL subtype sequence).
- Performance (reported result: copies/ml, IU/ml, pg/ml, ranges, %CD4 of lymphocytes, single/dual platform).
- Laboratory infrastructure requirements (high/low throughput, footprint, skill/training).
- GCLP (good clinical laboratory practise)
- Applies to equipment, assay, reagent and even sample collection/handling.





A typical approach to evaluation

- **Phase I (background and set up)**
 - Select appropriate comparative technology (Gold standard, more than one assay, more than one site, use automation)
- **Phase II (Design and analysis)**
 - Sample size: A balance of cost, risk of taking too few samples to measure lack of agreement.
 - Include reference material and controls
 - Statistical analysis (continuous data – can convert into discrete/bin approach):
intra/inter variability=background variability)





A protocol for method comparison

- Describe and summarize the data
- Visualize the data
- Choose the correct model for method comparison



- Analyze the data in sequence of sample preparation





- **Phase III (Reporting)**

- Validation report
- Good documentation
- Store everything!!!
- Take action if deviations

- **Phase IV (Follow-up)**

- Handle change control
- Participate in
EQA/proficiency testing

Clinically acceptable differences

- **CD4:**

- ~20cells/ul @ 200cells/ul (NB: data range)

- **Viral load:**

- 0.3 log copies/ml for intra-variability
- 0.5 log copies/ml for inter-variability
- 1.0 log copies/ml = clinical difference/patient mismanagement

Experiences

■ CD4:

- TetraCHROME (Beckman Coulter)
- *PanLeucogated* CD4 (Beckman Coulter)
- Flow Count (Beckman Coulter)
- Easy CD4 (Guava Technologies)
- FACSCount (Becton Dickinson)
- PointCare/AuRICA (PointCare Technologies)

■ Viral load

- COBAS (Roche) (Ampliprep/Amplicor)
- TaqMan (Roche)
- EasyQ (bioMerieux) (miniMAG/easyMAG)
- LUX assay (WITS, in house)
- P24 (Perkin Elmer)
- RT (Cavidi)

lscott@witshealth.co.za

Soon to launch methodcomparison.com

